

# TECHNICAL DATA SHEET DRAPE PERCUTANEOUS

(EU) 2017/745 Annex XI-Part A Production Quality Assurance.

# **Drape PERCUTANEOUS Features**

| Product Description         | Protective equipment used by doctors or nurses for surgical operations to prevent the possible risk of infection.   |  |  |  |  |
|-----------------------------|---|--|--|--|--|
| Product Class               | (EU) 2017/745 Medical Device Regulation – Class Business Rule I   |  |  |  |  |
| Manufacturer'<br>s Location | Tio Medikal 2/20 st. No.:53 (Begos 3. North<br>Entrance, 35400 Buca OSB/Buca/İzmir  |  |  |  |  |
| Purpose of usage            | The product is used for kidney operations. There is a 25x25cm window in the middle of the drape and it is covered with a polyurethane incision film. There is a fluid collection bag under the incision |  |  |  |  |

## Quality

• Manufactured under ISO 13485:2016 and 13795-1 quality management standards.

area.

• It has CE certificate.

## **Bio-Compatibility**

- Does not contain latex.
- Sterilized with ethylene oxide.

#### **Related Standard**

- ISO 13485:2016 Medical Medical Devices Quality Management System
- 9001:2015 Quality Management System
- TS EN 13795 Surgical Clothing and Drapes
- TS EN ISO 11607-1:2017 Packaging for finally sterilised medical devices Part 1: Rules for materials, sterile barrier systems and packaging systems / Products are made in accordance with the relevant standard.

• TS EN ISO 11607-2:2017 Packaging of medical devices-Finally sterilised-Part 2: Validation requirements for forming, sealing and joining processes/ Products are made in accordance with the relevant standard.

### **Shelf Life**

• 3 years

## 1. Drape Percutaneous Sizes




Figure 1. Drape percutaneous Technical Demonstration

### **Available Materials:**

01 SMS (Spunbond PP/Meltblown/Spunbond PP) ,35-43 gr/m2 02 Dublex(PE /Viscose or PE/Spunbond PP) , 54-56 gr/m2

DIMENSIONS (CM)

**NUMBER IN COLUMN** 

REF. CODES

A

В

50×80×50

40×60×40

422.0X.001.01

300

150


75

30

**Tolerances:** -+3% (ALL SIZES ARE PRODUCED & CUSTOMISATION)

## 2. Percutaneous Pack Content

| Percutaneous Pack        | Size      | QTY |
|--------------------------|-----------|-----|
| Percutaneous drape       | 150*300cm | 1x  |
| Back table cover         | 150*200cm | 1x  |
| Op tape                  | 10*30*cm  | 1x  |
| Fluoroscopy cover        |           | 1x  |
| Camera cover             | 14*250cm  | 1x  |
| Reinforced surgical gown | L         | 2x  |
| Medical towel            | 40*40cm   | 2x  |
| Sterilisation drape      | 100*100cm | 1x  |



Page 3 / 5
Addex Medical

|   |                | EN 13795 requirements   |              |           |   |
|---|----------------|---|--------------|-----------|---|
| Characteristic                                  | Test Method    | High performance  Critical product Less critical product area |              | Results   | Main fabric   |
| Resistance to microbial penetration - Dry (CFU) | EN ISO 22612   | Not required  | <300         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Resistance to microbial penetration - Wet (/ B) | EN ISO 22610   | 6,0   | Not required | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Cleanliness - Microbial<br>(CFU/100cm²)         | EN ISO 11737-1 | <300  | <300         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Cleanliness - Particulate<br>matter (IPM)       | EN ISO 9073-10 | <3,5  | <3,5         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Linting (Log10 (lint count))                    | EN ISO 9073-10 | <4,0  | <4,0         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Resistance to liquid penetration                | EN 13795       | ≥ 100   | ≥ 10         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Bursting strength - dry                         | EN 13795       | ≥ 40  | ≥ 40         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Bursting strength - wet                         | EN 13795       | Not required  | ≥ 40         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Tensile strength - dry                          | EN 13795       | ≥ 20  | ≥ 20         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |
| Tensile strength - wet                          | EN 13795       | Not required  | ≥ 20         | Conformed | - SS Hydrophobic Spunbond<br>- SS Hydrophylic Spunbond<br>- Surgigal drape and gowns-Blue/white non woven<br>- Surgigal drape and gowns-Blue laminated non<br>woven |

## **BIOCOMPATIBILITY**

Biocompatability studies has been done, and test results are given in the below table. Test results are within limits.

Table 31 Biocomtability Study Documentation

| Test Name       | Test Method  | Test facility | Report No | Report<br>Date | Result      |  |
|-----------------|--------------|---------------|-----------|----------------|-------------|--|
| Cytotoxicity    | ISO10993-    | HACETTEPE     | ARGEDS-   | 01.08.2016     | Confirmed   |  |
|                 | 5:2010       | UNIVERSITY    | 2016/37   | 01.08.2016     | Confirmed   |  |
| Sensitizasyon   | ISO10993-10: | HACETTEPE     | ARGEDS-   | 24.10.2016     | Confirmed   |  |
| Sensitization   | 2014         | UNIVERSITY    | 2016/37   | 24.10.2016     | Confirmed   |  |
| Cilt İrritasyon | ISO10993-10: | HACETTEPE     | ARGEDS-   | 11.06.2016     | Confinenced |  |
| Skin Irritation | 2014         | UNIVERSITY    | 2016/37   | 11.06.2016     | Confirmed   |  |

#### Instruction for use:

- 1. The package of product shall be opened in sterile and aceptic conditions.
- 2. For a clean peel, open the package from the direction of the arrow slowly.
- 3. There are labels or marks on the drapes, which helps the user to follow the patient direction and unfolding directions.
- 4. For the is incision film or adhesive tape on the drapes, first of all peel the paper carrier and fix the drape on the operational area of the patient. After then unfold the drape by following the directions.
- 5. Surgical drapes are ready for the operation, when they are completely unfolded.

| <b>®</b> | Do not use if package is damaged       |     | Manufacturer          | <b>(2)</b> | Single Use<br>(Do not re-use)  | Ω              | Use by date               |
|----------|--|-----|-----------------------|------------|--|----------------|---------------------------|
| 类        | Do not expose the product to sunlight. | REF | Catalogue<br>number   | []i        | Consult instructions for use   | 30 °C          | Temperature<br>Limitation |
| STERINZE | Do not re-sterilize                    | ~   | Manufacturing<br>Date | STERILEEO  | Sterilized using<br>Ethyleneoxide and<br>Single sterile barier<br>system | <b>CE</b> 2696 | CE Marking                |
| *        | Keep Dry                               | LOT | Batch Code            | $\bigcirc$ | Caution  | MD             | Medical<br>device         |
| <u></u>  | Unique device identifier               |     |                       |            |  |                |                           |